




Guidance Input and Output (GIO) Module



Hardware Introduction and Reference Manual

Version 4.0.0, June 18, 2015

FINAL RELEASE

P/N: G1X0-DI-01102

Document Content

The information contained herein is the property of Precise Automation Inc., and may not be copied, photocopied, reproduced, translated, or converted to any electronic or machine-readable form in whole or in part without the prior written approval of Precise Automation Inc. The information herein is subject to change without notice and should not be construed as a commitment by Precise Automation Inc. This information is periodically reviewed and revised. Precise Automation Inc. assumes no responsibility for any errors or omissions in this document.

Copyright © 2004-2015 by Precise Automation Inc. All rights reserved.

The Precise Logo is a registered trademark of Precise Automation Inc.

Trademarks

GIO, GSB, Guidance 3400, Guidance 3300, Guidance 3200, Guidance 2600, Guidance 2400, Guidance 2300, Guidance 2200, Guidance 1400, Guidance 1300, Guidance 1200, Guidance 0200 Slave Amplifier, Guidance 0006, Guidance 0004, Guidance Controller, Guidance Development Environment, GDE, Guidance Development Suite, GDS, Guidance Dispense, Guidance Input and Output Module, Guidance Programming Language, GPL, Guidance Slave Board, Guidance System, Guidance System D4/D6, PrecisePlace 1300, PrecisePlace 1400, PrecisePlace 2300, PrecisePlace 2400, PreciseFlex 400, PreciseFlex 1300, PreciseFlex 1400, PrecisePower 300, PrecisePower 500, PrecisePower 2000, PreciseVision, RIO are either registered or trademarks of Precise Automation Inc., and may be registered in the United States or in other jurisdictions including internationally. Other product names, logos, designs, titles, words or phrases mentioned within this publication may be trademarks, service marks, or trade names of Precise Automation Inc. or other entities and may be registered in certain jurisdictions including internationally.

Any trademarks from other companies used in this publication are the property of those respective companies. In particular, Visual Basic, Visual Basic 6 and Visual Basic.NET are trademarks of Microsoft Inc.

Disclaimer

PRECISE AUTOMATION INC., MAKES NO WARRANTIES, EITHER EXPRESSLY OR IMPLIED, REGARDING THE DESCRIBED PRODUCTS, THEIR MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS EXCLUSION OF IMPLIED WARRANTIES MAY NOT APPLY TO YOU. PLEASE SEE YOUR SALES AGREEMENT FOR YOUR SPECIFIC WARRANTY TERMS.

Precise Automation Inc.
727 Filip Road
Los Altos, California 94024
U.S.A.
www.preciseautomation.com

Warning Labels

The following warning and caution labels are utilized throughout this manual to convey critical information required for the safe and proper operation of the hardware and software. It is extremely important that all such labels are carefully read and complied with in full to prevent personal injury and damage to the equipment.

There are four levels of special alert notation used in this manual. In descending order of importance, they are:



DANGER: This indicates an imminently hazardous situation, which, if not avoided, will result in death or serious injury.



WARNING: This indicates a potentially hazardous situation, which, if not avoided, could result in serious injury or major damage to the equipment.



CAUTION: This indicates a situation, which, if not avoided, could result in minor injury or damage to the equipment.

NOTE: This provides supplementary information, emphasizes a point or procedure, or gives a tip for easier operation

Table Of Contents

Introduction to the Hardware	1
Guidance Input and Output (GIO) Module Overview	1
Installation, Configuration and Software	2
Mounting and Installation	2
GIO Hardware and Software Configuration	2
Accessing GIO Remote DIO	5
Hardware Reference	6
Guidance Input and Output (GIO) Hardware Overview	6
Digital Input and Output Signals	7
RS-485 Signal / 24VDC Power Connector	10
RS-485 Termination Jumper	11
Status Red/Green LED	12
Unit Number Jumpers	12
GIO Board Jumpers	13
Appendix A: Product Specifications	16
Guidance Input and Output (GIO) Module Specifications	16

Introduction to the Hardware

Guidance Input and Output (GIO) Module Overview

The Guidance Input and Output (GIO) module is a very compact device that connects to any Guidance Controller and adds an additional 12 optically isolated remote digital input channels and 8 optically isolated remote digital output channels.

GIO modules interface to a controller via a two-wire, bi-directional, daisy chained RS-485 line and can be located up to approximately 6 meters from the controller. Depending upon the timing requirements of the application and the available 24VDC power, a mix of as many as 8 GIO's or Guidance Slave Boards (GSB's) may be connected to a RS-485 cable. (For PreciseFlex 400 Sample Handlers with Linear Rails that include two internal GSB's, only 2 additional GIO's or GSB's can be connected due to 24VDC power limitations.)



WARNING: The GIO contains unshielded 24 VDC signals and pins. This product is intended to be mounted in a cabinet or machine chassis that is not accessible when power is turned on.


Once the GIO and the Guidance Controller are interfaced, the on/off states of the GIO digital inputs and outputs are automatically transferred to the controller via the Precise Servo Network protocol. Application programs running on the controller can access the GIO input and output values as though they were local to the controller. The only consideration is that GIO values are delayed somewhat since they are periodically transferred via the network protocol.

For systems that require an even greater number of remote digital inputs and/or outputs or other types of IO interfaces or that require that the I/O be located a greater distance from a Guidance Controller, Precise sells an Ethernet-based Remote Input and Output Module (RIO). For detailed information on the Guidance Controller or the RIO, please see the *Precise Documentation Library*.

Installation, Configuration and Software

Mounting and Installation

The Guidance Input and Output module is an open frame device with exposed 24VDC and must be mounted inside of a cabinet or other enclosure. The mounting holes are shown in **Blue** in the following drawing with dimensions in millimeters. There should be sufficient airflow across this board to ensure that the components do not become excessively hot.



WARNING: The GIO contains unshielded 24 VDC signals and pins. This product is intended to be mounted in a cabinet or machine chassis that is not accessible when power is turned on.

GIO Hardware and Software Configuration

GIO Board Unit Number

Up to a maximum of 8 GIO's or Guidance Slave Boards (GSB's) can be theoretically interfaced to a Guidance Controller. Therefore, this section describes how to configure and address up to 8 GIO or GSB

boards. However, due to communication timing and 24VDC power considerations, which are described below, **a maximum combination of 4 GIO's or GSB's is a more practical limit except for special system configurations.**

GIO's and GSB's can be connected to the RS-485 daisy chain in any order without altering their operation or identification. A GIO or GSB is not identified by its physical position in the daisy chain but by a unit number that is embedded in its low-level communication messages. This unit number is specified by jumpers on the GIO. A GIO's unit number is automatically combined with "GIO_" to generate a keyword that is used to configure the communication protocol between the GIO and its master controller. The Unit number can be arbitrary selected and need not be sequentially assigned, but each GIO or GSB must have a unique unit number within a given RS-485 system.

The setting of the GIO's [Unit Number Jumpers](#) is presented in the following table along with the corresponding unit number and keyword identifier.

NOTE: The Unit Number Jumpers for the GIO are J7/J8/J9 whereas the GSB utilizes J8/J9/J10.

J7	J8	J9	GIO Unit	GIO Keyword
In	In	In	1	GIO_1
Out	In	In	2	GIO_2
In	Out	In	3	GIO_3
Out	Out	In	4	GIO_4
In	In	Out	5	GIO_5
Out	In	Out	6	GIO_6
In	Out	Out	7	GIO_7
Out	Out	Out	8	GIO_8

RS-485 Signal Termination

There is one hardware configuration option that is dependent upon the ordering of modules in the RS-485 daisy chain. For noise immunity, **termination jumpers must be installed on the GIO's or GSB's or controller on the extreme ends of the RS-485 daisy chain.** The termination jumpers must be removed for all controllers or boards in between. On the GIO board, the [Termination Jumper](#) is labeled J6. Consult the hardware description for your specific master controller to determine its RS-485 termination jumper location.

Controller Software Configuration

In order for the controller to communicate with an IO module, the GIO's GIO_Keyword must be entered into the "Servo network node identifier" (DataID 151) parameter database array in the master controller. This provides the controller with the information it requires to communicate with the GIO board. The position of the GIO_Keyword in the DataID 151 array assigns the board a "network node number". Within the controller's software environment, the network node number (and not the GIO's unit number) is used to reference DIO on the GIO board. By convention, the first network node is always the master controller and the first element of the DataID 151 is always the controller's serial number.

Guidance IO (GIO) Manual

For example, if the GIO Unit Number Jumpers are set to select unit #4 (Keyword GIO_4), to define the GIO as the second network node, DataID 151 should be defined as follows:

```
DataID 151: "<master>", "GIO_4", "", "", "", "", ...
```

Normally, only elements at the end of the DataID 151 array can be blank and a blank entry cannot precede any non-blank servo node entry. As a special case, blank entries may occur before a GIO keyword. So you could enter DataID 151 as:

```
DataID 151: "<master>", "", "", "", "", "GIO_4"
```

This allows GIO boards to be assigned to fixed node numbers regardless of how many other slave boards might be in the system. Having a fixed node number means the I/O signal numbers are also fixed.

24VDC Logic Power Considerations

The GIO must be provided with 24VDC to power both the board's logic and to drive the input and output signals. The amount of 24VDC power available for GIO boards may limit the number of boards that can be wired in a system.

As a wiring convenience, GIO's and GSB's typically draw the 24VDC power from the same 10-pin daisy chained ribbon cable that provides the RS-485 signals. In this configuration, the Guidance Controller and its associated 24VDC power supply provide the logic and signal power to the GIO's and GSB's.

Starting in early 2013, all Guidance Controllers can output a maximum of 2A at 24VDC on the 10-pin RS-485 connector assuming that the controller's 24VDC power supply has sufficient power. Prior to early 2013, this was limited to only 1.35A.

The minimum power requirement for the GIO's logic is 0.05A. In the worst case where a board's digital outputs are all configured as sourcing and are driving 100mA and the digital inputs are configured as sourcing as well, a single GIO could draw 1.05A or more. Even with 2A available from the controller, at most, two GIO's could be supported by the power available from a controller.

Fortunately, in a typical system, even if the GIO's digital outputs are configured as sourcing, 20mA to 50mA of drive per channel is more common. So, we would expect a typical GIO configured with sourcing outputs to draw 0.5A, which would permit up to 4 GIO's to be interfaced to a single controller.

If GSB's are interfaced as well, their power consumption must be taken into account as well as any motor power that is drawn from the controller's 24VDC via the RS-485 cable.

If the controller's available power on the RS-485 cable is not sufficient for an application, an external 24VDC power source can be wired to the GIO and GSB boards (although this is not as convenient) or the digital signals could be configured as sinking if practical.

Communication Timing Considerations

Each GSB exchanges messages with the master controller every motion control trajectory cycle. The trajectory period is determined by the parameter "Trajectory Generator update period in sec" (DataID 600) on the master controller. This parameter typically ranges from 1 to 4 msec and determines the delay in reading or writing input and output values on the GIO.

As the number of GIO and/or GSB boards increases, the trajectory period must also be increased in order to accommodate the increased transmission times. The table below shows the maximum number of GSB or GIO nodes possible for different trajectory periods.

Trajectory period (msec)	Maximum number of GIO or GSB nodes
1	1
2	4
4	8

Accessing GIO Remote DIO

Once the GIO and the Guidance controller are configured, the general-purpose digital input and output signals on the GIO can be easily accessed from the controller as remote I/O signals. The update period for these signal values is the same as the trajectory period set by "Trajectory Generator update period in sec" (DataID 600) on the master controller.

To specify a GIO's general-purpose digital I/O signal from the master controller, multiply the GIO's network node number (not the GIO's unit number) by 100000 and add the I/O signal's number. For example, to access a GIO board configured as the 2nd or nth network node, add 200000 or 100000*n to the signal number.

Node	Signal Offset	GIO Module	
		Outputs (8)	Inputs (12)
2	200000	200013-200020	210001-210012
n	100000*n	n00013-n00020	n10001-n10012

Hardware Reference

Guidance Input and Output (GIO) Hardware Overview


The Guidance Input and Output module provides additional digital input and output signal interfaces that further enhance a Guidance Controller's ability to connect to external equipment. The interfaces and configuration hardware for this module include the following:

- [Digital Inputs / Outputs Connector](#)
- [RS-485 Signal / 24VDC Power Connector](#)
- [RS-485 Termination Jumper](#)
- [Status Red/Green LED](#)
- [Unit Number Jumpers](#)

The picture below illustrates the top surface of the GIO and identifies each of the user connectors and the major configuration components. To jump to the detailed information for a specific connector, click on the connector interface name or the connector.




WARNING: The GIO contains unshielded 24 VDC signals and pins. This product is intended to be mounted in a cabinet or machine chassis that is not accessible when power is turned on.



In the following sections, the pin outs for each of the connectors and the settings for the jumpers are described.


Digital Input and Output Signals

The GIO board provides 12 general purpose optically isolated digital input signals and 8 general purpose optically isolated digital output signals. These signals are presented in a single 26-pin IDC compatible connector with retaining latches. This type of connector permits these signals to be easily interfaced to other devices.




Depending upon the setting of the [Sinking Versus Sourcing Inputs Jumpers \(J3, J4, J5\)](#), the inputs are configured as "sinking" or "sourcing" in groups of 4 signals. J3 configures inputs 1 through 4, J4 configures inputs 5 through 8, and J5 configures inputs 9 through 12.

If an input signal is configured as "sinking", the external equipment must provide a 5VDC to 24VDC voltage to indicate a logical high value or no voltage for a logical low. This configuration is compatible with "sourcing" (PNP) sensors.



If an input signal is configured as "sourcing", the external equipment must pull the signal input pin to ground to indicate a logical high and must let the line float high to 24VDC to signal a logical low value. This configuration is compatible with "sinking" (NPN) sensors.

SOURCING DIGITAL INPUT




As shipped from the factory, all digital inputs are normally configured as "sourcing".


By setting eight sets of [Sinking Versus Sourcing Output Jumpers \(J11\)](#), each output can be individually configured as either "sinking" or "sourcing".

If an output signal is "sinking", the external equipment must provide a 5VDC to 24VDC pull-up voltage on the output pin and the GIO pulls this pin to ground when the signal is asserted as true. This configuration is compatible with "sourcing" (PNP) devices.

SINKING DIGITAL OUTPUT



If an output signal is "sourcing", the external equipment must pull-down the output pin to ground and the GIO pulls this pin to 24VDC when the signal is asserted as true. This configuration is compatible with "sinking" (NPN) devices.


SOURCING DIGITAL OUTPUT

As shipped from the factory, all digital outputs are normally configured as "sinking".

The pin out for the Digital Input and Output Connector and the corresponding GPL signal numbers are described in the following table. ***For the GPL signal numbers "n" is the GIO's Network Node number that is defined when configuring the controller, not the GIO's unit number.***

Signal	GPL Signal Number	Pin	Pin	GPL Signal Number	Signal
GND		1	2		GND
Input 1	n10001	3	4	n10002	Input 2
Input 3	n10003	5	6	n10004	Input 4
Input 5	n10005	7	8	n10006	Input 6
Input 7	n10007	9	10	n10008	Input 8
Input 9	n10009	11	12	n10010	Input 10
Input 11	n10011	13	14	n10012	Input 12
24VDC		15	16		24VDC
Output 1	n00013	17	18	n00014	Output 2
Output 3	n00015	19	20	n00016	Output 4
Output 5	n00017	21	22	n00018	Output 6
Output 7	n00019	23	24	n00020	Output 8
24VDC		25	26		GND
User Plug Parts	Molex 90142-0026. The sockets for this plug are 90110-2110 and the Molex crimp tool is a 63811-1000. Alternately, an IDC plug can be used, Samtec part # HCSD-13-01-N				

If the GIO board is installed within an enclosure, it is often convenient to use a ribbon cable whose ends are terminated with an IDC plug and a female DB25 connector to present the digital input and output signals on the enclosure's facilities panel. For example, if a GIO is installed inside a PreciseFlex 400 Sample Handler, this type of cable is used to connect the GIO to a DB25 plug that can be accessed from the outside at the base of the robot.



In this case, the 26th pin of the GIO's Digital Input and Output Signals Connector is not utilized. Also, the pin numbers for the signals are changed due to the DB25 pin numbering scheme. The revised pin numbers and the corresponding GPL signal numbers are described in the following table. **For the GPL signal numbers "n" is the GIO's Network Node number that is defined when configuring the controller, not the GIO's unit number.**


Signal	GPL Signal Number	Pin	Pin	GPL Signal Number	Signal
GND		1	14		GND
Input 1	n10001	2	15	n10002	Input 2
Input 3	n10003	3	16	n10004	Input 4
Input 5	n10005	4	17	n10006	Input 6
Input 7	n10007	5	18	n10008	Input 8
Input 9	n10009	6	19	n10010	Input 10
Input 11	n10011	7	20	n10012	Input 12
24VDC		8	21		24VDC
Output 1	n00013	9	22	n00014	Output 2
Output 3	n00015	10	23	n00016	Output 4
Output 5	n00017	11	24	n00018	Output 6
Output 7	n00019	12	25	n00020	Output 8
24VDC		13			
User Plug Parts	DB25 Male Plug				

RS-485 Signal / 24VDC Power Connector

The GIO communicates with a Guidance Controller using a RS-485 interface. RS-485 is a two-wire, bi-directional, multi-drop, daisy chained, high-speed serial interface. Once the GIO and controller are connected and configured, the controller's operating system automatically manages reading and writing the digital input and output signals on the GIO at a regular interval without requiring special programming.

The RS-485 signals are provided in an IDC connector. To simplify wiring, this connector also provides the 24VDC and ground lines that are necessary to operate the GIO. If a single GIO module is interfaced

to a Guidance Controller, a simple ribbon cable with an IDC connector on each end can connect the GIO to the controller and provide both communication signals and power to the GIO.



For reliable communications, **termination jumpers must be installed on the GIO's, GSB's or controllers on the extreme ends of the RS-485 daisy chain.** The termination jumpers must be removed for all controllers or boards in between.

The pin out for the RS-485/Power Connector is described in the following table.

Pin	Description
1	24 VDC. A minimum of 0.05 Amps is required for the GIO's logic power. A maximum of 1 Amp additional is required when all 8 digital outputs are configured as sourcing and are driving 100mA each and all of the digital inputs are configured as sourcing as well. If the digital outputs are driving less than 100ma each, the additional 1 Amp will be reduced accordingly. In a typical system, sourcing outputs normally drive 20mA to 50mA.
2	
3	GND
4	GND
5	GND
6	GND
7	GND
8	RS485+
9	RS485-
10	GND
User Plug Part No	AMP 746285-1 or Molex 22-55-2101 or 90142-0010. For the Molex plug, use Molex pins 16-02-0103 and Molex crimp tool 63811-1000.

RS-485 Termination Jumper

In order for the RS-485 daisy chained serial bus to operate properly, the ends of the bus must be electrically terminated. This electrical termination prevents transmitted signals from being reflected back into the cable and corrupting valid data. However, interior boards in the daisy chain must not have any electrical termination.

To allow a GIO to be placed anywhere in a RS-485 daisy chain, this board includes configurable bus termination that is controlled by [two jumper posts at position J6](#) on the top surface of the board.



When a GIO is placed at either end of a RS-485 daisy chain, the two posts must be jumpered together to terminate the bus. When a GIO is placed at an interior node of a chain, the jumper must be removed. **As shipped from the factory, the jumper is installed and the GIO is ready to be connected at either end of a chain.**

Status Red/Green LED

The GIO module includes a Status LED mounted on its top surface. This indicates the power and execution state of the board. The interpretation of the red and green LED is described in the following table.

LED State	System Status	Description
Continuously off or on	(1) Logic power off or (2) CPU crashed	Normally indicates that 24VDC logic power is off. In rare instances, indicates that the GIO CPU has crashed due to a system hardware or software error.
Alternating red and green	Board booting	The 24VDC logic power has been turned on and the GIO board is being initialized. <i>If this state continues for more than a minute or two, it usually indicates a hardware failure.</i>
Blinking red	Board operating but not communicating	The GIO CPU has completed its startup process and is operational. However, RS-485 communication with the master Guidance Controller has not been established.
Blinking green	Normal operation, IO active	The GIO is operational and is actively communicating with the master Guidance Controller. The input and output signals on the board are active.

Unit Number Jumpers

In the low-level RS-485 communications, the "unit number" determines which GIO is the originator or recipient of each message, not the position of the GIO board in the RS-485 daisy chain. This unit number is configurable using a group of three jumper posts on the top of the GIO board. The unit numbers can be arbitrarily assigned and do not have to be sequential, but they do have to be unique within a controller system.



The unit number also determines a keyword ("GIO_<unit_number>") that is specified to configure a GIO board as a node in a controller's Servo Network. **Note, at the software application level, the network node number and not the GIO board unit number determines [how the GIO's digital input and output lines are addressed](#).**

For more information on node numbers and configuring the controller, please see the [Hardware and Software Configuration](#) section of this manual.

In the following table, the interpretation of the [Unit Number Jumpers \(J7, J8, J9\)](#) is provided. **As shipped from the factory, all of the jumpers are installed, which indicates GIO unit #1.**

NOTE: The Unit Number Jumpers for the GIO are J7/J8/J9 whereas the GSB utilizes J8/J9/J10.


J7	J8	J9	GIO Unit	GIO Keyword
In	In	In	1	GIO_1
Out	In	In	2	GIO_2
In	Out	In	3	GIO_3
Out	Out	In	4	GIO_4
In	In	Out	5	GIO_5
Out	In	Out	6	GIO_6
In	Out	Out	7	GIO_7
Out	Out	Out	8	GIO_8

GIO Board Jumpers

The GIO board has a number of hardware jumpers that determine the configuration of various hardware functions. Depending upon the type of jumper, there may be two, three or five jumper posts. Posts are tied (shorted) together using black jumper plugs. The five wide jumper posts for configuring the digital output signals are shown below.



The locations of each of the key sets of jumpers are illustrated in the following diagram and are identified by stenciled labels on the surface of the GIO board.



The following table describes each of the sets of jumpers and how the pins must be shorted ("jumped") in order to set a specific configuration. When a direction (e.g. left versus right) is described, it is with respect to the GIO board oriented as shown in the picture above.

Jumpers	Description	Setting
J11 Digital Outputs Sink/Source	These jumpers determine if each of the Digital Output Signals is "sinking" or "sourcing". This array of jumpers consists of eight columns of five posts. Each column determines the setting for a single digital output. The left column (pins 1-5) control to the first digital output signal. To set an output to sourcing, two jumpers must be installed on the top four posts. To set an output to sinking, two jumpers must be installed on the bottom four posts. For example, to set the first output to sinking, a jumper should short pins 2 & 3 and a second jumper should short pins 4 & 5. As shipped	For sinking, short bottom four posts using two jumpers For Sourcing, short top four posts using two jumpers

	from the factory, all outputs are set to sinking.	
J10 Spare Jumper	This is the right most jumper in the J7/J8/J9/J10 group and is currently unused. As shipped from the factory, this jumper always is installed.	Always installed
J7/J8/J9 Unit Number	In the low-level RS-485 communication protocol, the Unit Number determines which GIO is the originator or recipient of each message, not the position of the GIO board in the RS-485 daisy chain. See the Unit Number Jumper section for a description of these jumpers. The left most jumper is J7. As shipped from the factory, all three of these jumpers are installed and the board is set to unit #1.	Install or remove to define GIO Unit Number.
J6 RS-485 Bus Termination	This jumper controls if RS-485 Bus Termination is enabled on this board. For reliable communications, if a GIO is at the end of a RS-485 daisy chain, this jumper must be installed to terminate the communication line. If a GIO is in the middle of a RS-485 daisy chain, this jumper must be uninstalled to disable the termination. As shipped from the factory, this jumper is installed and the GIO should be installed at the end of the RS-485 daisy chain.	Install jumper J6 to terminate the RS-485 communication lines.
J3/J4/J5 Digital Inputs Sink/Source	These jumpers configure the "sinking" or "sourcing" mode of operation for the Digital Input Signals in groups of four. J3 configures inputs 1 through 4, J4 configures inputs 5 through 8, and J5 configures inputs 9 through 12. Each jumper consists of three posts. If the left posts are jumpered, the group of input signals is configured as sinking. If the right posts are jumpered, the group of inputs is configured as sourcing. As shipped from the factory, all inputs are set to sourcing.	For Sinking, jumper left two posts. For Sourcing, jumper right two posts.
Status LED	This is a green and red LED that blinks to indicate the operational status of the controller.	

Appendix A: Product Specifications

Guidance Input and Output (GIO) Module Specifications

General Specification	Range & Features
Interface to Guidance Controller	
Communications Interface	Interfaces via a two-wire, bi-directional, daisy chained RS-485 line and can be located up to approximately 6 meters away from the controller
Communications Protocol	Operates as part of the Precise Servo Network.
Scanning Rate	Input and output states are updated at the rate set by the "Trajectory Generator update period in sec" (DataID 600) of the master controller. This update rate is typically 1-4 msec.
Number of units	A combination of up to 8 GIO's and GSB's can theoretically be simultaneously interfaced to a Guidance Controller. The actual maximum is a function of the "Trajectory Generator update period in sec" (DataID 600) of the master controller and the available 24 VDC power. In typical systems, a maximum of 4 GIO or GSB boards can be simultaneously operated.
Input and Output Interfaces	
Digital Input Channels	12 optically isolated digital inputs Configurable in banks of 4 as sinking or sourcing 5VDC to 24VDC for logic high if sinking 24VDC supplied for logic high if sourcing
Digital Output Channels	8 optically isolated digital outputs Individually configured as sinking or sourcing 24VDC maximum pull up
General	
Size and Weight	50mm (W) x 70mm (L) x 15.5mm (H), 0.040 kg
Low Voltage Logic Power	24VDC required for logic and input/output functions
	A minimum of 0.05 Amps is required for logic power. A maximum of 1 Amp additional is required when all 8 digital outputs are configured as sourcing and are driving 100mA each and all of the digital inputs are configured as sourcing. If the digital outputs are driving less than 100ma each, the additional 1 Amp will be reduced accordingly. In a typical system, sourcing outputs normally drive 20mA to 50mA.