

## Details:

96 Well Glass Bottom Microplate | 0.17mm | Low Glass  
 630µl | Black Plastic Upper | Standard Height Clear Lid

### Dimension Specification (mm)

A .....	14.38
B .....	11.24
C .....	127.76
D .....	113.38
E .....	85.48
F .....	74.24
Overall Height .....	14.35
Glass Height .....	0.88
Well to Well Spacing .....	9.00
Well Depth .....	11.50
Glass Thickness .....	0.17
Well Volume .....	630µl

